User manual Heart rotation 27apr2022


1. Enter the rotation list to rotate and/or view current status

Enter the 'Heart rotation' menu through the donor tab.




When you enter the 'heart rotation' system you will find:

- The current status




- The place to make an actual rotation




2. Rotation

- 1. Choose which center you wish to rotate (=receiving center)
- 2. Select cause of rotation
- 3. *Optional* enter a comment e.g. donor Scandia number. (You do not have to enter the order of the center as this is automatically logged and is found in the rotation log).



When you have registered the correct information click on the 'Rotate' button



Now the displayed rotation list will be update